



Instructions for the *AeroChamber MAX*[®] Valved Holding Chamber (VHC)

AeroChamber
MAX[®]

Below, please find abbreviated instructions, but please remember to read the complete instructions that are packaged with your **AeroChamber MAX**[®] VHC device.



- 1 Remove cap from the MDI and chamber. *Before use, ensure the instructions supplied with the MDI and chamber have been read.*
- 2 Shake the MDI immediately before each use. Insert the MDI into the backpiece of the chamber.
- 3 Put mouthpiece into mouth and close lips around it to ensure an effective seal. The **Flow-Vu**[®] Inspiratory Flow Indicator only moves if the patient has a good seal.
- 4 Breathe out gently and depress the MDI at the beginning of a slow inhalation. Use the **Flow-Vu**[®] Inspiratory Flow Indicator to assist in the coordination of this maneuver.
- 5 Breathe in slowly and deeply through the mouth until a full breath has been taken.
- 6 Hold breath for 5-10 seconds, if possible. Otherwise keep lips tight on the mouthpiece to maintain seal for 2-3 breaths after MDI is depressed. **Administer one (1) puff at a time.**
- 7 Slow down inhalation if you hear the **FLOWSignal**[®] Whistle sound.
- 8 Follow instructions supplied with MDI on how long to wait before repeating instructions as prescribed.



monaghan[™]